

# Mobile REMM for Windows Mobile - README

Version 1.5

This README will be available in the Mobile REMM group of your PC's start menu when the installation completes.

## **Contents**

- System Requirements
- Installation
  - Quick-Start Instructions
  - Detailed Instructions
  - Upgrade Instructions
  - Resolving Installation Errors
- Uninstalling
- Known Issues
- Release History
- Mobile REMM Overview
- Disclaimer

## **System Requirements**

To install and run Mobile REMM you need the following:

- A Windows Mobile "Pocket PC" (touch screen) device or Smartphone.
- Pocket PC 2003 (version 4.x) or later devices are supported (includes Windows Mobile 2003, Windows Mobiles 5.0 – 6.5)
- Approximately 3.0MB of storage space, either on the device or on an expansion card, for Mobile REMM
- Internet Explorer mobile (formerly Pocket Internet Explorer) must be present on the device.
- A personal computer with the following installed software:
  - Device synchronization software. On Windows XP, Microsoft ActiveSync 4 or later is required. On Windows Vista & Windows 7, the required Windows Mobile Device Center is part of Vista & Windows 7, but it must be upgraded to version 6 or later. For the latest versions of either, go to <http://www.microsoft.com/windowsmobile/activesync/default.mspx>.

To determine the version of the device's operating system, select Start→Settings, select the System tab, and then select About.

To determine the amount of memory on your device, from the same System tab, select Memory, and then the Main tab. If your device has an expansion card, you can go to the Storage Card tab to determine its free space.

## ***Installation***

If a previous version of Mobile REMM is already installed, see the "*Upgrade Instructions*" section below.

The Mobile REMM installation consists of two parts: installation to your desktop/laptop PC (documentation and files required by ActiveSync), and installation to your Windows Mobile device. After the PC portion of the installation, your device synchronization software will download Mobile REMM to your Windows Mobile device and launch its installation. If your device is not in its cradle, then the Mobile REMM installation will continue the next time your device is cradled.

## **Quick-Start Instructions**

1. Cradle your device.
2. Download the setup program (MREMM\_WinMobile\_Setup.exe) to your PC and run it.
3. The PC files, including this README will be installed.
4. The Windows Mobile setup program will start and install Mobile REMM for Windows Mobile onto your device.
5. Select the default location on the device, or select an alternative location (e.g., a storage card).
6. Upon completion of the installation, run Mobile REMM by selecting the Mobile REMM icon at Start→Programs.

## **Detailed Instructions**

1. Download to your PC the MREMM\_PocketPC\_Setup.msi setup program from the REMM website and run it (double-click it). This installs the Mobile REMM documentation to your PC, as well as the files that Microsoft ActiveSync requires for installation to the device. By default, these files are located under the folder "C:\Program Files\Mobile REMM\PocketPC", and consist of:

```
C:\Program Files
    \Mobile REMM
        \PocketPC
            README.rtf
            Disclaimer.txt
        \DeviceInstallerFiles (files downloaded to the device during
synchronization)
```

2. When the desktop installation completes, the Windows Mobile installation is launched. This will install Mobile REMM for Windows Mobile application.

The remaining details of the installation depend on whether you have a Windows Mobile 5.0 – 6.5 or Pocket PC 2003 device.

Regardless, upon completion of the installation, run Mobile REMM by selecting the Mobile REMM icon at Start→Programs.

## **Windows Mobile 5.0 – 6.5**

After completion of the PC portion of the installation, respond to any prompts presented on you

Windows Mobile device:

1. If your device has a storage card, you will be prompted on the device to choose a location (Device or Card) to install Mobile REMM for Windows Mobile – note that optimum performance will be realized by installing to the device
2. An unknown publisher warning for the “Mobile REMM for Windows Mobile.CAB” may be presented. Answer “Yes”.
3. You may be prompted to restart the device.

Mobile REMM is installed under the folder “\Program Files\Mobile REMM” on either the device or storage card, depending on what was selected for the destination.

### **Pocket PC 2003**

1. A dialog will appear on your PC with the message *“Install Mobile REMM and other pending applications using the default application install directory?”*. By default, the Mobile REMM application files will be installed in your device’s main memory at “\Program Files\Mobile REMM”. If space on your device is limited, you may need to consider placing Mobile REMM on an expansion card instead, if your device is so equipped. Note that optimum performance will be realized by installing to main memory.
  - a. Answer “Yes” to continue installation to the default location if your device meets the memory requirements, OR
  - b. Answer “No”. A “Select Destination Media” dialog is presented with a list of options. The default option will be “Main Memory”, and an indication of the amount of required and available space will be provided. If your device has an expansion card, it will be included as an option in this list, typically as “Storage Card”. Select the desired option and then press OK. On a storage card, Mobile REMM is installed under a folder “\Mobile REMM for Windows Mobile”
2. The Windows Mobile installer will proceed to download Mobile REMM to your device.
3. At completion of the download, an “Application Downloading Complete” dialog will be displayed on your PC and present the following message: “Please check your device for any additional operations that must be performed on the device to complete installation.” This is an indication that your device may be presenting additional dialogs to complete the installation.

### **Upgrade Instructions**

An upgrade is performed the same as a normal installation, but with the following considerations.

1. Windows Mobile 5.0 – 6.5 devices: For Mobile REMM for Windows Mobile, a dialog will be present indicating that the previous version will be removed before the new one is installed. Select “OK”.
2. Pocket PC 2003 devices: If there is not enough space on the device, this will be indicated on the Select Destination Media dialog by the Space Available and Space Required values. Press OK, and a “Not Enough Space on Device dialog” will be presented. Select OK to indicate that you want to continue anyway. The calculation of available space does not, unfortunately, account for the fact that the previous version of Mobile REMM will be removed, and thus that space is available. Therefore, continuation of the installation should be okay.
3. As with a normal install, if the device displays prompts asking if files should be overwritten, answer “Yes to All”.

## Resolving Installation Errors

None.

## *Uninstalling*

### Desktop Uninstall

The desktop portion of the Mobile REMM installation can be uninstalled from your PC via the Control Panel's Add/Remove Programs option.

### Windows Mobile Uninstall

There are two ways to uninstall Mobile REMM from your Windows Mobile:

1. On the device, go to Start→Settings→System→Remove Programs. Select Mobile REMM from the list of programs, and then select Remove.
2. If your device is cradled with an ActiveSync connection, open ActiveSync and select Tools→Add/Remove Programs. A dialog listing the applications on your device is displayed. Find Mobile REMM and clear the checkbox.

## *Known Issues*

None.

## *Release History*

Application V1.5.4:

- Updated algorithm for diagnosing and managing radiation contamination: how to manage bioassay specimens
- New section: Understanding Radiation Units

Application V1.4.0

- Updated the look and feel

Application V1.3.0:

- Changed Banner.

Application V1.2.0:

- Updated the Management Algorithm
- Updated the Dose Estimator

Application v1.1.0:

- Updated Dose Estimator

Application v1.0.0:

- Initial public release

## *REMM Overview*

- REMM is a tool that helps in diagnosis and treatment of victims of mass casualty radiological/nuclear events
- Mobile REMM has only selected pages from REMM Web site - <http://remm.nlm.gov>
- See REMM Web site for much more data, images, and background material
- See REMM Web site for key "printable pages" to carry in your pocket

- See REMM downloadable version for access without the internet

Mobile REMM currently exists as a standalone application for

- Windows Mobile PDAs (Pocket PC and Smartphone)
- iPhone, iPod and iPad
- Android devices
- BlackBerry Smartphone.

### ***Feedback***

Contact us at [nlmremm@mail.nih.gov](mailto:nlmremm@mail.nih.gov) with questions/feedback.

### ***Disclaimer***

The U.S. Government does not warrant or assume any legal liability or responsibility for the accuracy, completeness, or usefulness of any information, apparatus, product, or process disclosed.

It is not the intention of NLM to provide specific medical advice to the public, but rather to provide users with information to better understand their health.

NLM does not endorse or recommend any commercial products, process, or services.